


Start your micro lateral flow application today



Multiplex lateral flow immunoassays - our complete solution

You can find everything you need for LFA assays in our offering, including buffers, membranes and more

NON-CONTACT DISPENSING


Dot or line printing of capture probes on LFA membranes



With our instruments or as a service inhouse

CONJUGATION

Detector molecules & spraying on conjugate pad



ASSAY INCUBATION

Biomarker detection using μ LFA cartridges



ARRAY IMAGING & ANALYSIS

sciREADER LF2 for qualitative/quantitative analysis of biomarkers



Our service packages to success - no matter your starting point

Provide us with your idea, sample, capture and detection molecules


MINIATURIZATION



Printing of dots on prestructured membranes

1 week


MULTIPLEXING



Combination of your singleplex LFA tests

3 weeks

PLATFORM TRANSFER



Conversion of microplate ELISA to μ LFA

3 weeks

FEASIBILITY OF YOUR OWN IDEA




Conceptional study for biomarker detection

4 weeks



CONTACT US



scienion.com