SCIENION
BY BICO

cellenONE® BSC

PRECISION & PROTECTION FOR AUTOMATED SINGLE CELL CLONING & OMICS.



About cellenONE ®

cellenONE® is a unique platform combining high accuracy single cell isolation and precision nanoliter reagent dispensing, allowing miniaturization of a range of single cell omics protocols and enabling automated cell line development workflows with best-in-class clonal outgrowth.

Benefits

- Integrated Class II biosafety protection
- Reliable monoclonality with imaging
- Multi-parameter fluorescence & brightfield sorting
- Sorts cells from 0.5–80 µm
- All-in-one single-cell sorting & isolation
- Nanoliter dispensing for cost-efficient workflows
- Hundreds of cells isolated in minutes
- Compatible with 96-, 384-, and 1536-well plates


Gene Editing

Stem Cell Reprogramming

Synthetic Biology

Monoclonal Antiboody Discovery

Microbiology



CONTACT US

SCIENION GmbH
Wagner-Régeny-Str. 15
12489 Berlin, Germany
Fon +49 (0)30 6392 1700
Fax +49 (0)30 6392 1701
support@scienion.com
www.scienion.com

SCIENION US, Inc
1231 W. Warner Road, Ste. 103
Tempe, AZ 85284
United States
Tel: +1 (888) 988-3842
USsalessupport@scienion.com

SCIENION (UK) Ltd
2000, Lakeside North Harbour
Western Road, Portsmouth
PO6 3EN
United Kingdom
+44 (0)7483 388 271
+44 (0)23 9323 3603
support@scienion.com

SCIENION China
Room 1106B Building 3
No.391 Guiping Road Xuhui
District, Shanghai
Phone: +86 186 2170 9292
support@scienion.com
www.scienion.com